

Can I use an in-house standard instead of a Ph. Eur. CRS?



You are here:

[EDQM FAQs](#) / [EDQM FAQs in English](#) / [PRODUCTS AND SERVICES \(Publications, Reference standards, CombiStats\)](#) / [EDQM reference standards](#) / [General information on EDQM reference standards](#) / Can I use an in-house standard instead of a Ph. Eur. CRS?

Answer:

In-house standards may be used for routine analysis, provided that they are calibrated against EDQM reference standards. For more information, see Ph. Eur. general chapter 5.12. *Reference standards*.