

How can I find out the purity of a Ph. Eur. CRS?



You are here:

[EDQM FAQs](#) / [EDQM FAQs in English](#) / [PRODUCTS AND SERVICES \(Publications, Reference standards, CombiStats\)](#) / [EDQM reference standards](#) / [Establishment of chemical reference standards \(CRS\)](#) / How can I find out the purity of a Ph. Eur. CRS?

Answer:

When required by the prescribed test and/or assay of the corresponding Ph. Eur. text, the content assigned a CRS is given in the leaflet available in the [Ph. Eur. reference standards database](#).

For impurity CRSs used for quantification in the related substances test of the corresponding Ph. Eur. monograph, the content can be assumed to be 100%, unless otherwise indicated in the leaflet.

In all other cases, no value is assigned and no assumption can be made on the purity or content of the CRS.